

Sling & Harness Guide tollos

Single Patient Use

Combi Slings

- Seated sling with head support & padded legs
- Comfortable soft poly blend fabrics
- Designed for use from the bed, chair, or floor
- Single patient use for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T015CPP-PSS	Combi Sling, Extra Small	600 lbs/ 272 kg	10/ box
T015CP-PSS	Combi Sling, Small	600 lbs/ 272 kg	10/ box
T015C-PSS	Combi Sling, Medium	600 lbs/ 272 kg	10/ box
T015CX-PSS	Combi Sling, Large	600 lbs/ 272 kg	10/ box
T015CXX-PSS	Combi Sling, XL	1000 lbs/ 454 kg	10/ box
T015CXXX-PSS	Combi Sling, XXL	1000 lbs/ 454 kg	10/ box
T015CXXXX-PSS	Combi Sling, XXXL	1000 lbs/ 454 kg	5/ box

 $Note: The\ combi\ sling\ is\ not\ intended\ for\ above-the-knee\ amputees.\ Tollos\ recommends\ the\ Long\ Seat\ sling\ instead.$

Long Seat Slings

- Seated sling designed for general lifting
- Comfortable soft poly blend fabrics
- Ideal for amputees or groin injuries
- Single patient use for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T015LPP-PSS	Long Seat Sling, Extra Small	600 lbs/ 272 kg	10/ box
T015LP-PSS	Long Seat Sling, Small	600 lbs/ 272 kg	10/ box
T015L-PSS	Long Seat Sling, Medium	600 lbs/ 272 kg	10/ box
T015LX-PSS	Long Seat Sling, Large	600 lbs/ 272 kg	10/ box
T015LXX-PSS	Long Seat Sling, XL	1000 lbs/ 454 kg	10/ box
T015LXXX-PSS	Long Seat Sling, XXL	1000 lbs/ 454 kg	10/ box
T015LXXXX-PSS	Long Seat Sling, XXXL	1000 lbs/ 454 kg	5/ box





Single Patient Use

Repositioning Slings

- Full-length, reclining sling
- Comfortable soft poly blend fabrics
- Designed for repositioning or turning in bed
- Perfect for floor rescues, even with a backboard in place
- Single patient use for infection control

Part No.	Description	Safe Working Load	Box Quantity
T017LX-PSS	Repo Sling, Standard	600 lbs/ 272 kg	10/ box
T017LXB-PSS	Repo Sling, Bariatric	1000 lbs/ 454 kg	10/ box

Multi-Purpose Bands

- Comfortable padded bands
- Soft poly blend fabrics
- Allow virtually clear access for wound care & hygiene procedures
- Single Patient Use for infection control
- Box includes 10 bands, 2 carabiners, & 2 straps

Part No.	Description	Safe Working Load	Box Quantity
T019R-PSS T019RO-PSS	Multi-purpose bands, 5.5" (140 mm) width 5.5" bands with no carabiners/ straps	200 lbs/ 91 kg	10/ box
T019RB-PSS	Multi-purpose bands, 12" (305 mm) width	400 lbs/ 181 kg	10/ box

Steady Aid® Harnesses

- Padded harness
- Comfortable soft poly blend fabrics
- Designed for sit-to-stand lifting & training
- Single patient use for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T016DP-PSS	Steady Aid [®] Harness, Small	350 lbs/ 158 kg	10/ box
T016D-PSS	Steady Aid [®] Harness, Medium	350 lbs/ 158 kg	10/ box
T016DX-PSS	Steady Aid [®] Harness, Large	450 lbs/ 204 kg	10/ box
T016DXX-PSS	Steady Aid [®] Harness, XL	450 lbs/ 204 kg	10/ box
T016DXXX-PSS	Steady Aid [®] Harness, XXL	750 lbs/ 340 kg	10/ box





Combi Slings

- Seated sling with head support & padded legs
- Comfortable, yet durable, poly blend fabrics or mesh
- Designed for use from the bed, chair, or floor
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T015CPPF T015CPPMF	Combi Sling, Extra Small Combi Sling, Mesh, Extra Small	600 lbs/ 272 kg	1 each
T015CPF T015CPMF	Combi Sling, Small Combi Sling, Mesh, Small	600 lbs/ 272 kg	1 each
T015CF T015CMF	Combi Sling, Medium Combi Sling, Mesh, Medium	600 lbs/ 272 kg	1 each
T015CXF T015CXMF	Combi Sling, Large Combi Sling, Mesh, Large	600 lbs/ 272 kg	1 each
T015CXXF T015CXXMF	Combi Sling, XL Combi Sling, Mesh, XL	1000 lbs/ 454 kg	1 each
T015CXXXF T015CXXXMF	Combi Sling, XXL Combi Sling, Mesh, XXL	1000 lbs/ 454 kg	1 each
T015CXXXXF T015CXXXXMF	Combi Sling, XXXL Combi Sling, Mesh, XXXL	1000 lbs/ 454 kg	1 each

Note: The combi sling is not intended for above-the-knee amputees. Tollos recommends the Long Seat sling instead.

Hygiene Elite Slings

- Seated sling designed for toiletting & peri care
- Comfortable, yet durable, poly blend fabrics
- Built-in extra back support for added comfort
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T015HPPF	Hygiene Elite Sling, Extra Small	600 lbs/ 272 kg	1 each
T015HPF	Hygiene Elite Sling, Small	600 lbs/ 272 kg	1 each
T015HF	Hygiene Elite Sling, Medium	600 lbs/ 272 kg	1 each
T015HXF	Hygiene Elite Sling, Large	600 lbs/ 272 kg	1 each
T015HXXF	Hygiene Elite Sling, XL	1000 lbs/ 454 kg	1 each
T015HXXXF	Hygiene Elite Sling, XXL	1000 lbs/ 454 kg	1 each



Long Seat Slings

- Seated sling designed for general lifting
- Comfortable, yet durable, poly blend fabrics or mesh
- Ideal for amputees or groin injuries
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T015LPPF T015LPPMF	Long Seat Sling, Extra Small Long Seat Sling, Mesh, Extra Small	600 lbs/ 272 kg	1 each
T015LPF T015LPMF	Long Seat Sling, Small Long Seat Sling, Mesh, Small	600 lbs/ 272 kg	1 each
T015LF T015LMF	Long Seat Sling, Medium Long Seat Sling, Mesh, Medium	600 lbs/ 272 kg	1 each
T015LXF T015LXMF	Long Seat Sling, Large Long Seat Sling, Mesh, Large	600 lbs/ 272 kg	1 each
T015LXXF T015LXXMF	Long Seat Sling, XL Long Seat Sling, Mesh, XL	1000 lbs/ 454 kg	1 each
T015LXXXF T015LXXXMF	Long Seat Sling, XXL Long Seat Sling, Mesh, XXL	1000 lbs/ 454 kg	1 each
T015LXXXXF T015LXXXXMF	Long Seat Sling, XXXL Long Seat Sling, Mesh, XXXL	1000 lbs/ 454 kg	1 each

Repositioning Slings

- Full-length, reclining sling
- Comfortable soft poly blend fabrics or mesh
- Designed for repositioning or turning in bed
- Perfect for floor rescues, even with a backboard in place
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Qty
T017LX T017LXM	Repo Sling, Large Repo Sling, Mesh, Large	600 lbs/ 272 kg	1 each
T017LXMG	Repo Sling, Mesh, Large, w/ elastic corner straps	600 lbs/ 272 kg	1 each
T017LXB T017LXMB	Repo Sling, XL Repo Sling, Mesh, XL	1000 lbs/ 454 kg	1 each
T017LXMBG	Repo Sling, Mesh, XL w/ elastic corner straps	1000 lbs/ 454 kg	1 each
T017LXBE	Repo Sling, XXL	1000 lbs/ 454 kg	1 each



100	00 lbs/ 454 kg	1 each
100	00 lbs/ 454 kg	1 each
100	00 lbs/ 454 kg	1 each
	Safe Working Load	Box Qty
	600 lbs/ 272 kg	1 each
ans	600 lbs/ 272 kg	1 each

Steady Aid® Harnesses

- Padded harness
- Comfortable soft poly blend fabrics
- Designed for sit-to-stand lifting & training
- Washable for infection control
- Color-coded binding



Note: Use of this product, or any sling product, should be done with the oversight of a trained medical professional. Do not use these slings without a clinician present.

Steady Aid Elite Harnesses

- Padded harness
- Comfortable soft poly blend fabrics
- Designed for sit-to-stand lifting & training
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
T016PDE T016PDEB	Steady Aid [®] Elite Harness, Small Harness with Gluteal Strap, Small	350 lbs/ 158 kg	1 each
T016DE T016DEB	Steady Aid [®] Elite Harness, Medium Harness with Gluteal Strap, Med	350 lbs/ 158 kg	1 each
T016DEL T016DELB	Steady Aid [®] Elite Harness, Large Harness with Gluteal Strap, Large	450 lbs/ 204 kg	1 each
T016DELX T016DELXB	Steady Aid [®] Elite Harness, XL Harness with Gluteal Strap, XL	450 lbs/ 204 kg	1 each
T016DELXX T016DELXXB	Steady Aid [®] Elite Harness, XXL Harness with Gluteal Strap, XXL	750 lbs/ 340 kg	1 each



Note: Use of this product, or any sling product, should be done with the oversight of a trained medical professional. Do not use these slings without a clinician present.

Walking Harnesses

- · Walking harness designed for gait training
- Helps protect from fall or injury
- Soft poly blend fabrics
- Washable for infection control
- Color-coded binding

Part No.	Description	Safe Working Load	Box Quantity
TSC015WSP	Walking Harness w/ Leg Straps, Small	350 lbs/ 158 kg	1 each
TSC015WS	Walking Harness w/ Leg Straps, Medium	350 lbs/ 158 kg	1 each
TSC015WSX	Walking Harness w/ Leg Straps, Large	450 lbs/ 204 kg	1 each
TSC015WSXX	Walking Harness w/ Leg Straps, XL	450 lbs/ 204 kg	1 each
TSC015WSXXX	Walking Harness w/ Leg Straps, XXL	750 lbs/ 304 kg	1 each

Multi-Purpose Bands

- Comfortable padded bands
- Soft poly blend fabrics
- Allow virtually clear access for wound care & hygiene procedures
- Washable for infection control
- Pair includes 2 bands, 2 carabiners, & 2 straps

Part No.	Description	Safe Working Load	Box Quantity
T019R T019R4	Set of two bands, 5.5" (140mm) width 5.5" bands with 4 carabiners/ 4 straps	200 lbs/ 91 kg	1 pair
T019RB T019RB4	Set of two bands, 12" (305mm) width 12" bands with 4 carabiners/ 4 straps	400 lbs/ 181 kg	1 pair

Dual Turner

- Comfortable soft poly blend fabrics
- Designed for turning patients in bed
- Washable for infection control

Part No.	Description	Safe Working Load	Box Quantity
T019DT T019DTM	Dual Turner, Standard Dual Turner, Mesh, Standard	350 lbs/ 158 kg	1 each
T019DTX T019DTXM	Dual Turner, Large Dual Turner, Mesh, Large	600 lbs/ 272 kg	1 each





General Information

Cleaning & Disinfecting

Reusable

Spot Cleaning: Most disinfection solutions are suitable for spot cleaning the reusable slings & harnesses. Bleach products can be used in a dilution of 100ppm or 0.01%, however, greater concentrations may cause premature deterioration of the fabrics.

Reusable slings should be laundered between patients.

Washing Instructions:

Recommended wash cycle: BS 5651:1978 Hospital Laundry Wash (Normal)

Recommended detergent: ECE Reference Detergent + Sodium Perborate + TAED Recommended Wash Cycle: 6 gallons (24 liters) of water at 165°F (74°C) for 15 minutes

Temp not to exceed 185°F (85°C).

Recommended Rinse Cycle: 4 cycles at 104°F (40°C) for 3 minutes per cycle Recommended Dry Cycle: Drying cycle: BS EN ISO 6330:2012 Procedure E

Low Temp, not to exceed 120°F (50°C)

Bleach products may be used during the wash cycle in a dilution of 100ppm or 0.01%.



Single-Patient-Use

Do not wash: single patient use. Discard sling after use.

Tollos single-patient-use sling labels include a wash indicator. If the sling has been washed and recycled for use again, the wash indicator will turn red. Do not use a sling with a red wash indicator.



Single-patient-use slings & harnesses may be spot cleaned with disinfecting solutions. Tollos does not recommend the use of bleach in any dilution on single-patient-use slings & harnesses.

Small 75 - 125 lbs Red 600 lbs/ 272 kg Medium 125 - 175 lbs Yellow 600 lbs/ 272 kg	Sizing Guide - F	Safe Working Load		
Medium 125 - 175 lbs Yellow 600 lbs/ 272 kg	X-Small	49 - 75 lbs	White	600 lbs/ 272 kg
	Small	75 - 125 lbs	Red	600 lbs/ 272 kg
174 050 W 0	Medium	125 - 175 lbs	Yellow	600 lbs/ 272 kg
Large 1/4 - 250 lbs Green 600 lbs/ 2/2 kg	Large	174 - 250 lbs	Green	600 lbs/ 272 kg
X-Large 250 - 400 lbs Blue 1000 lbs/ 454 kg	X-Large	250 - 400 lbs	Blue	1000 lbs/ 454 kg
XX-Large 400 - 600 lbs Orange 1000 lbs/ 454 kg	XX-Large	400 - 600 lbs	Orange	1000 lbs/ 454 kg
XXX-Large 600 - 1000 lbs Gray 1000 lbs/ 454 kg	XXX-Large	600 - 1000 lbs	Gray	1000 lbs/ 454 kg

Transfer Aids

Sally Tubes & Sheets

- Assists in moving patients easily
- Helps protect clinicians and providers from injury
- Washable reusables or single-patient-use for infection control



Part No.	Description	Box Quantity
ST-5829S	Sally Tube, Single Patient Use, 58 in x 29 in (145 cm x 73 cm) 10 each per box, 5 boxes per case	50/ case
ST-5848S	Sally Tube, Single Patient Use, 58 in x 48 in (145 cm x 125 cm) 10 each per box, 5 boxes per case	50/ case

www.tollos.com

Corporate Office 8 Easter Court, Suite J Owings Mills, MD 21117 410.363.1515 PHONE 888.363.7224 TOLL FREE

Central Zone Office 930 N Beltline Rd Suite 107 Irving, TX 75061 888.363.7224 TOLL FREE

Western Zone Office 5560 Skylane Blvd Suite C Santa Rosa, CA 95403 888.363.7224 TOLL FREE

Canadian Office 75 Dyment Road Barrie, ON L4N 3H6 705.733.0022 PHONE 800.565.7075 TOLL FREE

