



Your patients need to move. You need to move them without putting yourself or them at risk.

Now it's easy to help your patients and residents move around. The Steady Aid Sit-to-Stand helps caregivers perform smooth seat-to-seat transfers.

The Steady Aid is easy to use and safe for you and your patients. It is specifically designed for seat-to-seat transfers and makes them virtually effortless.

Steady Aid Sit-to-Stand

tollos

The Steady Aid® Sit-to-Stand is one of the industry's safest, most ergonomic patient lifts on the market. Manufactured in North America, this lift includes a long list of standard features specifically engineered to optimize patient comfort and safety, as well as reduce caregiver injury.

The Steady Aid offers a smooth, fluid transition from a sitting position to a standing position without a lot of fuss.

Steady Aid Features

- Safe working load: 500 lb
- Smart Technology
- Power base widening
- Removable foot plate
- EZ Roll casters
- Integrated head protector & knee pad
- Four-function hand controller
- Two harness attachment points

Standard equipment includes battery operation, durable construction, emergency up/down switch, emergency stop switch, and powered base-widening system, allowing for integration with virtually every bed, chair, or stretcher.

Supplement your lift with Tollos' Steady Aid Elite Harnesses.

TECHNICAL SPECIFICATIONS

Safe Working Loads:	Standard - 500lbs (227 kg)
Operation:	Electric Actuator
Power Source:	Two 12 V, 2.6 Ah batteries
Charger:	24 V CSA/UL approved
Construction:	Rectangular Steel Tubing
Finish:	Polyester powder coat, baked
Overall Length:	35 inches (889 mm)
Overall Weight:	132 lb (60 kg)
Boom - Highest Position:	60 inches (1524 mm)
Boom - Lowest Position:	22 inches (559 mm)
Base Width - Closed:	25 inches (635 mm)
Base Width - Open:	38 inches (965 mm)
Front Casters:	3 inches (76 mm)
Rear Locking Casters:	4 inches (102 mm)
Knee Pad:	Multi-position, adjustable
Emergency Stop Switch:	Standard
Emergency Up/Down Switch:	Standard



The Steady Aid bariatric sit-to-stand offers more safety and security for your bariatric patients.



The Steady Aid can be your extra hands in helping your patients up to a standing position.



The Steady Aid Elite Harness offers added support to your patient while transferring seat-to-seat.

