



Lifting and moving patients or residents has always been a challenge.

With the Ultralift Mobile Lift, you can lift, move, reposition, and turn your patients without strain to you or fear of injury for them.

Protect yourself and protect those you care for.

Ultralift offers an easy, powerful option and it is ideally suited for turning, gait training, rehabilitation, repositioning, and lifting, even from the floor.

---

# Ultralift Mobile Lift

tollos

The Ultralift® Mobile Lift is one of the industry's safest, most user-friendly patient lifts on the market. Manufactured in North America, this lift includes a long list of standard features specifically engineered to optimize patient comfort and safety, as well as reduce caregiver injury. The Ultralift is also suitable for effective gait training and general rehabilitation activities.

## Ultralift Mobile Lift Features

- Smart Technology
- Power base widening
- Extra long boom
- EZ Roll casters
- True six-point spreader bar
- Extended lifting range of 56 inches
- Optional beavertail legs available

Standard equipment includes a powered base-widening system, battery operation, durable construction, emergency stop switch, emergency up/down switch, and our true, six-point spreader bar.



Optional beavertail legs fit under low beds and gurneys, while still allowing for safe, comfortable transport.



An ergonomic control panel with a dual battery warning system and dual controls means the Ultralift is always charged and ready to use.

## TECHNICAL SPECIFICATIONS

Safe Working Load:	750 lbs (340 kg)
Operation:	Electric Actuator
Power Source:	Two 12 V, 2.6 Ah batteries
Charger:	24 VCSA/UL approved
Construction:	Body - Rectangular Steel Tubing Legs - Rectangular Steel Tubing
Finish:	Polyester powder coat, baked
Overall Weight:	235 lb (107 kg)
Height (top of mast):	52.5 inches (1334 mm)
Spreader Bar (Highest):	78.6 inches (1996 mm)
Spreader Bar (Lowest):	22 inches (559 mm)
Base Width - Closed:	24.5 inches (622 mm)
Base Width - Opened:	42 inches (1067 mm)
Total Leg Length:	48.75 inches (1238 mm)
Working Leg Length:	34.25 inches (870 mm)
Turn Radius:	54.25 inches (1378 mm)
Front Casters:	3 inches (76 mm) (optional: 4", 102 mm or low-profile beavertail 2", 51mm)
Rear Locking Casters:	4 inches (102 mm) (optional: 5", 127 mm or low-profile beavertail 3", 76mm)
Emergency Stop Switch:	Standard
Emergency Up/Down Switch:	Standard
External Scale:	Optional

Supplement your lift with Tollos' complete line of Slings and Harnesses.



The true six-point spreader bar with spring-loaded locking clips make this lift multi-purpose, including repositioning, rehabilitation, limb fitting, and turning.



Lift patients of up to 750 lbs, rotate a patient for an unobstructed 360 degrees, accomplish 'head first' floor lifts, and more with the Ultralift Mobile Lift.